

# ANATOMY

## LOWER EXTREMITIES



# LOWER EXTREMITIES

## NATURAL SIZE




**EI-28D**  
**Complete bony structure of the foot.**

It is supplied unmounted, the bones are loose.  
Dimensions: 25 x 10 x 6 cm.  
Sh. wt. 0.19 Kgs.



**EI-30D**  
**Complete bony structure of the foot, mounted.**

Threaded on wire.  
Dimensions: 25 x 10 x 6 cm. Sh. wt. 0.19 Kgs.



**EI-35D**  
**Complete bony structure of the lower extremity.**

Threaded on wire.  
Dimensions: 100 x 25 x 12 cm.  
Sh. wt. 0.60 Kgs.



**EI-36D**  
**Complete bony structure of the lower extremity. It reproduces half pelvis.**

Threaded on wire.  
Dimensions: 100 x 25 x 12 cm. Sh. wt. 0.90 Kgs.



**EI-31D**  
**Complete bony structure of the foot with tibia and fibula section.**

Threaded on wire.  
Dimensions: 25 x 10 x 15 cm. Sh. wt. 0.29 Kgs.

YOU CAN CHOOSE RIGHT OR LEFT EXTREMITIES, ADD "D" FOR RIGHT, "I" FOR LEFT.

# LOWER EXTREMITIES

## NATURAL SIZE



**ED-27D**  
**Complete bony structure of the foot, didactic.**

It is supplied unmounted, the bones are loose.  
Dimensions: 25 x 10 x 6 cm.  
Sh. wt. 0.19 Kgs.


**EP-29D**  
**Complete bony structure of the foot, mounted.**

Threaded on nylon, (not illustrated).  
Dimensions: 25 x 10 x 6 cm.  
Sh. wt. 0.19 Kgs.



**ER-33D**  
**Complete bony structure of the foot.**

Threaded on wire.  
Foot painted muscle origins and insertions  
Dimensions: 25 x 10 x 6 cm. Sh. wt. 0.19 Kgs.



**ER-34D**  
**Complete bony structure of the lower extremity.**

Threaded on wire.  
Lower extremity painted muscle origins and insertions.  
Dimensions: 100 x 25 x 12 cm.  
Sh. wt. 0.60 Kgs.

**ED-25D**  
**Complete bony structure of the foot, didactic, mounted.**

Threaded on nylon, (not illustrated).  
Dimensions: 25 x 10 x 6 cm. Sh. wt. 0.19 Kgs.

**ED-26D**  
**Complete bony structure of the foot, didactic, mounted.**

Threaded on wire, (not illustrated).  
Dimensions: 25 x 10 x 6 cm. Sh. wt. 0.19 Kgs.

**ER-39D**  
**Complete bony structure of the lower extremity. It reproduces half pelvis.**

Threaded on wire.  
Lower extremity painted muscle origins and insertions.  
Dimensions: 100 x 25 x 12 cm. Sh. wt. 0.90 Kgs.

YOU CAN CHOOSE RIGHT OR LEFT EXTREMITIES, ADD "D" FOR RIGHT, "I" FOR LEFT.

**Avda. Sistema Solar, 13 – Nave 11  
28830 San Fernando de Henares  
Madrid**

**Telf. +34 918 328 318  
Telf. móvil: +34 657 783 403**

**email: [info@simulacionessanitarias.com](mailto:info@simulacionessanitarias.com)  
<http://www.simulacionessanitarias.com>**